# ImAnOx

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•  
*Colorimetric test system for the determination of antioxidative capacity in serum and EDTA-plasma*

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---

<b>1. INTENDED USE</b>	<b>12</b>
<b>2. SUMMARY AND EXPLANATION OF THE TEST</b>	<b>12</b>
<b>3. PRINCIPLE OF THE TEST</b>	<b>12</b>
<b>4. MATERIAL SUPPLIED</b>	<b>13</b>
<b>5. MATERIAL REQUIRED BUT NOT SUPPLIED</b>	<b>13</b>
<b>6. PREPARATION AND STORAGE OF REAGENTS</b>	<b>14</b>
RECONSTITUTION OF REACTION BUFFER B	14
PREPARATION OF REAGENT 1	14
PREPARATION OF REAGENT 2A AND 2B	14
RECONSTITUTION OF CALIBRATORS AND CONTROLS	14
<b>7. PRECAUTIONS</b>	<b>15</b>
<b>8. SPECIMEN COLLECTION AND PREPARATION</b>	<b>15</b>
<b>9. ASSAY PROCEDURE</b>	<b>16</b>
PROCEDURAL NOTES	16
SAMPLE AND STANDARD PREPARATION	16
PIPETTING SCHEME	17
<b>10. RESULTS</b>	<b>17</b>
CALCULATION	17
<b>11. LIMITATIONS</b>	<b>18</b>
<b>12. QUALITY CONTROL</b>	<b>18</b>
EXPECTED VALUES	18
CONTROLS	18
<b>13. PERFORMANCE CHARACTERISTICS</b>	<b>19</b>
PRECISION AND REPRODUCIBILITY	19
DETECTION LIMIT	19
<b>14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE</b>	<b>20</b>

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•

*Colorimetric test system for the determination of antioxidative capacity in serum and EDTA-plasma*

Gültig ab/valid from 25.11.2001

Artikelnummer/Catalogue no.:

KC 5200

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96 Tests/96 determinations

Lagerung/Storage:

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2-8 °C, calibrator at -20 °C

CE

## 1. INTENDED USE

The *Immundiagnostik* Assay is intended for the quantitative determination of antioxidants in EDTA-plasma and serum. This Assay is designed for in vitro diagnostic use only.

## 2. SUMMARY AND EXPLANATION OF THE TEST

The human body is constantly under attack from free radicals that occur as part of normal cell metabolism, and by exposure to environmental factors such as UV light, cigarette smoke, environmental pollutants and gamma radiation. The resulting „Reactive Oxygen Species“ (ROS) circulate freely in the body with access to all organs and tissues, which can have serious repercussions throughout the body. The body possesses a number of mechanisms both to control the production of ROS and to cope with free radicals in order to limit or repair the damage to tissues

Overproduction of ROS or insufficient defence mechanisms lead to a dangerous disbalance in the organism. Thereby several pathomechanisms implicated in over 100 human diseases, e.g. cardiovascular disease, cancer, diabetes mellitus, inflammatory disease, aging, etc., were induced.

Determination of the **antioxidative capacity** becomes of fundamental importance in medical diagnosis and in research. The **ImAnOx**-Assay is fast, reliable and easy to perform. The **total antioxidative capacity** is measured.

## 3. PRINCIPLE OF THE TEST

The determination of the antioxidative capacity is performed by the reaction of antioxidants in the sample with a defined amount of exogenously provided hydrogen peroxide ( $H_2O_2$ ). The antioxidants in the sample eliminate a certain amount of the provided hydrogen peroxide. The residual  $H_2O_2$  is determined colorimetrically by an enzymatic reaction which involves the conversion of TMB to a colored product.

After addition of a stop solution, the samples are measured at 450 nm in a microtiter plate reader. The quantification is performed by the delivered calibrator.

The difference between applied and measured peroxide concentration in a defined time period is proportional to the reactivity of the antioxidants of the sample (antioxidative capacity). Quantification is performed with the enclosed calibrator.

**Please note:** As the reaction speed of the distinct antioxidants in the sample is different, the measured concentrations of antioxidants are equivalent to the reactivity of the distinct antioxidants and not to their total amount in the sample. Therefore, we use in our test system hydrogen-peroxide equivalents as unit for the antioxidative capacity.

#### 4. MATERIAL SUPPLIED

<i>Catalog No</i>	<b>Kit Components</b>	<b>Quantity</b>
KC5200ka	Calibrator (lyoph. 250 µl)	3 vials
KC5200ko	Control 1 and 2 (lyoph. 250 µl)	1 vial each
KC5200mp	Microtitre plate	2 pieces
KC5200pl	Peroxide solution	0.25 ml
KC5200ra	Reaction buffer A	105 ml
KC5200rb	Reaction buffer B (lyophilised)	2 vials
KC5200lc	Solution C	1.5 ml
KC5200el	Enzyme solution	35 µl
KC5200sl	ELISA stop solution, ready to use	11 ml
KC5200rl	Reconstitution solution (aqua bidest.)	5 ml

#### 5. MATERIAL REQUIRED BUT NOT SUPPLIED

- Vortex mixer
- Precision pipettors calibrated to deliver 10-100 µl
- ELISA-reader
- Incubation chamber for 37 °C

## 6. PREPARATION AND STORAGE OF REAGENTS

### Reconstitution of reaction buffer B

- Resuspend the reaction buffer B with **600 µl** solution C and mix for 2 minutes on a vortex-mixer. The reaction buffer B is light sensitive and should be stored in the dark. It is stable for 3 month at 2-8 °C.

This is just an example. The volumes are sufficient for 50 determinations. For larger or less approaches, the amounts of buffer have to be adjusted to the desired volumes.

### Preparation of reagent 1

- **Reagent 1 has to be prepared directly before use:**

Add **5 ml** reaction buffer A and **10 µl** peroxide solution = **dilution 1**

Add **100 µl** of dilution 1 and **4.9 ml** of reaction buffer A = **reagent 1**

Please note: The reagent 1 is not stable and can not be stored.

### Preparation of reagent 2a and 2b

- **Reagent 2a and 2b have to be prepared directly before use:**

Because of the self-absorption of the sample and the reaction products produced by the added hydrogen peroxide which are both detectable at 450 nm, it is of great importance to measure the sample **with** and **without** addition of enzyme.

- **Reagent 2a:**

Add **5 ml** reaction buffer A + **100 µl** reaction buffer B + **5 µl** enzyme solution **Important:** The enzyme solution must be sealed with parafilm after use.

- **Reagent 2b:**

Add **5 ml** reaction buffer A + **100 µl** reaction buffer B

### Reconstitution of calibrators and controls

The calibrators and controls have to be reconstituted in **250 µl** reconstitution solution (aqua bidest.).

After 5 minutes they should be homogenized on a vortex mixer. Take aliquots and store them at -20°C (the reconstituted calibrator and controls are stable for at least 2 weeks at -20°C). Avoid repeated freeze-thaw cycles. The concentration might have minor changes from lot to lot.

\*Additional controls can be ordered separately. Catalog No.: KC5200ko

## 7. PRECAUTIONS

- For in vitro diagnostic use only.
- This product contains human source material which was tested and found to be non-reactive to HBsAg, anti-HIV-1/2, and anti-HCV. Since no method can offer complete assurance that hepatitis B virus, HIV-1/2, HVC or other infectious agents are absent, these reagents should be handled as if potentially infectious.
- The stop solution contains acid. It must be handled with care. It can cause burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spill should be wiped out immediately with copious quantities of water. Do not breath vapor and avoid inhalation.
- Reagents should not be used beyond the expiration date shown on kit label.

## 8. SPECIMEN COLLECTION AND PREPARATION

### EDTA-plasma and serum

- Venous fasting blood is suitable for this test system. The blood sample can be shipped at 2-8 °C in at least 24 hours. Serum and EDTA-plasma should be stored at -20 °C up to the measurement. They are stable at -20 °C for 1 week.
- Lipaemia and haemolysis interfere with the test system. Such samples should not be measured.
- Samples with visible amounts of precipitates should be centrifuged (5 min at 10000 g) prior to measurement, and the resulting supernatant is used in the test.

## 9. ASSAY PROCEDURE

### Procedural notes

- The quality control guidelines should be observed.
- Incubation time, incubation temperature and pipetting volumes of the different components are defined by the producer. Any variations of the test procedure, that are not coordinated with the producer, may influence the test results. Immundiagnostik can therefore not be held reliable for any damage resulting from this.
- The assay should always be performed due to the manual which is given in the kit.

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**10 µl** sample, calibrator and control are pipetted in duplicates in appropriate wells.

**100 µl** reagent 1 is added.

Incubate for **10 min** at **37°C**.

**100 µl** reagent 2a respectively reagent 2b is added.

Incubate for **5 min** at **room temperature**.

Add **50 µl** stop solution.

**Measurement:** is performed immediately after addition of the stop solution at 450 nm in the ELISA reader.

## Pipetting scheme

(example)

	with enzyme		without enzyme		with enzyme		without enzyme		with enzyme		without enzyme	
	1	2	3	4	5	6	7	8	9	10	11	12
A	P1	P1	P1	P1	P9	P9	P9	P9	P14	P14	P14	P14
B	P2	P2	P2	P2	C1	C1	C1	C1	P15	P15	P15	P15
C	P3	P3	P3	P3	P10	P10	P10	P10	P16	P16	P16	P16
D	P4	P4	P4	P4	Cal	Cal	Cal	Cal	P17	P17	P17	P17
E	P5	P5	P5	P5	P11	P11	P11	P11	P18	P18	P18	P18
F	P6	P6	P6	P6	C2	C2	C2	C2	P19	P19	P19	P19
G	P7	P7	P7	P7	P12	P12	P12	P12	P20	P20	P20	P20
H	P8	P8	P8	P8	P13	P13	P13	P13	P21	P21	P21	P21

Cal: Calibrator

C1: Control 1

C2: Control 2

Px: Patient sample

**10. RESULTS**

## Calculation

The difference of the measurement of the sample „with enzyme“ and „without enzyme“ is inverse proportional to the antioxidative capacity:

For evaluation subtract the OD-values of the measurement „without enzyme“ from the OD-values of the measurement „with enzyme“ to get the  $\Delta$ OD. The antioxidative capacity can be calculated with the following formula.

$$\text{Antiox. capacity of the sample } (\mu\text{mol/l}) = 392 - (392 - \text{Conc. Calibrator}) * \frac{\Delta\text{OD sample}}{\Delta\text{OD Calibrator}}$$

## 11. LIMITATIONS

Whole blood or heparin-plasma can not be used.

Strong haemolytic and lipaemic samples often show pathological concentrations. We recommend not to measure those samples.

## 12. QUALITY CONTROL

Expected values

Normal values (preliminary):

EDTA-plasma and serum:

< 280 $\mu\text{mol/l}$	low antioxidative capacity
280-320 $\mu\text{mol/l}$	middle antioxidative capacity
> 320 $\mu\text{mol/l}$	high antioxidative capacity

The mean value of an apparently healthy group of people were found at 305  $\mu\text{mol/l}$  (n=69).

We recommend, that each laboratory should establish its own normal range. The values mentioned above are preliminary data and can deviate from other results published.

### Controls

Control samples or serum pools should be analyzed with each run of calibrators and patient samples. Results generated from the analysis of the control samples should be evaluated for acceptability using appropriate statistical methods. In assays in which one or more of the quality control sample values lie outside the acceptable limits, the results for the patient sample may not be valid.

### 13. PERFORMANCE CHARACTERISTICS

#### Precision and reproducibility

**Intra-Assay CV:** 0.9 % (181  $\mu\text{mol/l}$ ) [n = 12]

2.3 % (217  $\mu\text{mol/l}$ ) [n = 12]

**Inter-Assay CV:** 1.63 % (185  $\mu\text{mol/l}$ ) [n = 12]

2.43 % (217  $\mu\text{mol/l}$ ) [n = 12]

#### Detection limit

130  $\mu\text{mol/l}$

## 14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- This assay was produced and put on the market according to the IVD guidelines of 98/79/EC.
- The test components contain organic solvents. Contact with skin or mucous membranes has to be avoided.
- All reagents in the test package are to be used for research only.
- The reagents should not be used after the date of expiry stated on the label.
- Single components with different lot numbers should not be mixed or exchanged.
- The guidelines for medical laboratories should be observed.
- Incubation time, incubation temperature and pipetting volumes of the different components are defined by the producer. Any variations of the test procedure that are not coordinated with the producer may influence the results of the test. Immundiagnostik AG can, therefore, not be held reliable for any damage resulting from this.

**Short Protocol for the determination of antioxidative capacity in serum/EDTA-plasma (see page 6 for more details)**

Allow reagents to reach room temperature.  
Prepare ELISA reagents 1 and 2a/b (see page 4)

Pipette **10 µl** calibrators, samples and controls into each well

Add **100 µl** of reagent 1

Incubate for **10 min** at 37°C

Add **100 µl** reagent 2a, respectively reagent 2b

Incubate for **5 min** at room temperature

Add **50 µl** stop solution into each well

**Immediately measure the extinction at 450**