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PMN-Elastase ELISA Kit

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
PMN-Elastase ELISA KIT

*For the in vitro determination of PMN-Elastase in serum, plasma,
seminal plasma and stool*

Valid from 03.01.2008

K 6840

1. INTENDED USE

The described Enzyme-Linked-Immuno-Sorbent-Assay (ELISA) is intended for the quantitative determination of **PMN-Elastase** in serum, seminal plasma and stool. It is for *in vitro* diagnostic use only.

2. INTRODUCTION

PMN-Elastase from human polymorphnuclear granulocytes is a glycoprotein of 30 kDa which belongs to the group of serine proteases. Active **PMN-Elastase** is released from azurophil granula of neutrophil granulocytes after irritation or disintegration. The determination of the **PMN-Elastase** in stool is used to record inflammatory reactions in which neutrophils are involved. Especially in Crohn's disease the inflammatory process goes hand in hand with an increased phagocytic activity and the biological decay of the phagocytic cells, thus, leading to an increased release of **PMN-Elastase** and other lysosomal enzymes.

Indication

- Activation marker for Morbus Crohn
- Chronic joint inflammation
- Bacterial infection, sepsis

3. MATERIAL SUPPLIED

Catalogue No	Content	Kit Components	Quantity
K 6840MTP	PLATE	Microtiter plate, precoated	12 x 8 wells
K 6840WB	WASHBUF	ELISA wash buffer concentrate 10x	2 x 100 ml
K 6840EP	EXBUF	Extraction buffer concentrate 2.5 x	2 x 100 ml
K 6840A2	AB	Detection antibody (Second antibody, mouse anti-PMN-Elastase, monoclonal), concentrate	600 µl
K6840K	CONJ	Peroxidase-labeled antibody (goat-anti-mouse-POD)	15 ml
K 6840ST	STD	Calibrator, lyophilized (see specification for amount)	4 vials
K 6840KO	CTRL	Control, lyophilized (see specification for range)	4 vials
K 6840TMB	SUB	TMB substrate (Tetramethylbenzidine), ready to use	15 ml
K 6840AC	STOP	ELISA stop solution, ready to use	15 ml

4. MATERIAL REQUIRED BUT NOT SUPPLIED

- Bidistilled water (aqua bidest.)
- Laboratory balance
- Precision pipettors calibrated to deliver 10-1000 µl
- Covering foil for the microtiter plate
- Horizontal microtiter plate shaker
- Multi-channel dispenser or repeating dispenser
- Centrifuge capable of 3000 x g
- Vortex-Mixer
- Standard laboratory glass or plastic vials, cups, etc.
- Microtiter plate reader 450 nm

5. PREPARATION AND STORAGE OF REAGENTS

- To run the assay more than once, ensure that reagents are stored at conditions stated on the label. **Prepare only the appropriate amount necessary for each assay** (e.g. diluted antibody is not stable). The kit can be used up to 4 times within the expiry date stated on the label.
- Reagent with volumes **less than 100 µl** must be centrifuged before use to avoid loss of volume.
- The **ELISA wash buffer concentrate** (WASHBUF) must be diluted with aqua bidest. **1:10** before use (100 ml WASHBUF + 900 ml aqua bidest.). Crystals could occur due to high salt concentration in the stock solutions. The crystals must be redissolved at 37°C using a water bath **before dilution of the buffer solutions**. The buffer concentrate (WASHBUF) is stable at 2-8°C until the expiry date stated on the label. Diluted solutions can be stored at 2-8°C for 1 month.
- The **EXBUF** (extraction buffer concentrate) should be diluted with aqua bidest. **1:2.5** before use (90 ml concentrate + 135 ml aqua bidest.). Crystals could occur due to high salt concentration in the stock solutions. The crystals must be redissolved in a water bath at 37°C before dilution. The **EXBUF** (extraction buffer concentrate) is stable at 2-8°C until the expiry date stated on the label. Diluted buffer solution can be stored in a closed flask at 2-8°C for fore months.
- The **second antibody** (detection antibody, AB) must be diluted **1:20** in ELISA wash buffer (add 500 µl AB to 9.5 ml of ELISA wash buffer). The undiluted antibody (AB) is stable at **2-8 °C** until the expiry date given on the label. **Diluted antibody solution is not stable and can not be stored.**
- The **standard** (STD) must be reconstituted with **0.5 ml** of aqua bidest. (please see the product specification).

For the **PMN-Elastase** calibration curve dilute the concentrate (S1) as described in the following example:

S 1

S 2 = 200 µl S1 + 400 µl ELISA wash buffer

S 3 = 200 µl S2 + 400 µl ELISA wash buffer

S 4 = 200 µl S3 + 400 µl ELISA wash buffer

(See specification for range)

As zero calibrator use the ELISA wash buffer.

- The **control** (CTRL) must be reconstituted with **0,5 ml** of aqua bidest. (please see the product specification). Mix thoroughly by gentle inversion to insure complete reconstitution.

Reconstituted standards and control are not stable and can not be stored.

- All other test reagents are ready for use. The test reagents are stable up to the date of expiry (see label of test package) when stored at 2-8°C.

6. SAMPLE PREPARATION

Extraction step of the stool sample

Please use the **extraction buffer** as a sample extraction buffer.

1. We recommend the use of a stool sample preparation kit for dosing of 100 mg of stool sample (e.g. Sample preparation kit from Roche Diagnostics, Mannheim, Germany; cat # 745804). Add 5 ml buffer to the stool sample in the preparation kit.

Constant buffer volume: 5 ml

Constant dilution factor: 1 : 50

2. Alternatively, stool samples can be manually weighted within the range of 80 – 120 mg. Please note the exact sample amount!

- a. Add **5 ml** buffer to the stool sample independent of the sample amount.

Constant buffer volume: 5 ml

Variable dilution factor

The **dilution factor varies** depending on the sample amount which must be considered in the subsequent calculations. Please see the list of the correction factors below:

Weight [mg]	Buffer Volume [ml]
80	62.5
82	60.9
84	59.5
86	58.1
88	56.8
90	55.6
92	54.3
94	53.2
96	52.1
98	51.0
100	50

Weight [mg]	Buffer Volume [ml]
102	49.0
104	48.1
106	47.2
108	46.3
110	45.5
112	44.6
114	43.9
116	43.1
118	42.4
120	41.6

- b. The **buffer volume** for the individual samples **varies** depending on the sample amount (see table).

Variable buffer volume

Constant dilution factor: 1 : 50

Therefore, the same dilution factor can be used for all samples in the subsequent evaluation of the results .

Weight [mg]	Buffer Volume [ml]
80	4.0
82	4.1
84	4.2
86	4.3
88	4.4
90	4.5
92	4.6
94	4.7
96	4.8
98	4.9
100	5.0

Weight [mg]	Buffer Volume [ml]
102	5.1
104	5.2
106	5.3
108	5.4
110	5.5
112	5.6
114	5.7
116	5.8
118	5.9
120	6.0

Afterwards, mix the weighed stool sample with the buffer and vortex for 30 sec. Transfer ca. 1 ml stool suspension in an Eppendorf-cup and centrifuge for 5 min at 13000 rpm. For analysis, pipet **100 µl** of the supernatant per well. **The supernatant is not stable and can not be stored.** We recommend to weight fresh sample amount for a new assay, if the analysis should be repeated.

Seminal plasma

Seminal plasma should be stored at -20°C and defrosted immediately before use. Centrifuge the seminal-plasma samples for 5 minutes at 10000 rpm.

The samples should be diluted **1:10 - 1:20** in the ELISA wash buffer depending on the inflammatory status of the patient.

Serum and Plasma samples

Preanalytic handling

Significant differences in the PMN-Elastase levels can be observed due to different sample preparation procedures, e. g. up to 10-fold higher serum levels compared to the plasma PMN-Elastase concentrations. The reasons are as follows:

The granulocytes are activated during the serum clotting and release elastase granulocyte-activating markers. The time between serum collecting and analysis as well as repeated freeze-thaw cycles don't cause a PMN-Elastase concentration shift.

On the contrary, in the case of plasma samples, varying the time between sampling and analysis or the number of freeze-thaw cycles will cause variation in the observed PMN-Elastase levels. Therefore, the preanalytical conditions of plasma samples should be held constant. This is a general requirement independent of the used test-system.

Immundiagnostik recommends the use of serum samples for PMN-Elastase determinations.

Fresh collected blood should be centrifuged within one hour. If not assayed on the same day, it should be stored at -20 °C. Lipemic or hemolytic samples should be not analysed. Samples should be mixed well before assaying. We recommend to carry out duplicate analysis on each test sample.

Serum samples should be diluted **1:500** with the ELISA-wash buffer before assaying.

Plasma samples should be diluted **1:100** with the ELISA-wash buffer before assaying.

7. ASSAY PROCEDURE

Principle of the test

In a first incubation step, **PMN-Elastase** in the sample is bound to polyclonal sheep-anti-PMN-Elastase antibodies (in excess), which are immobilized on the surface of the microtiter wells (PLATE). To remove all unbound substances, a washing step is carried out. In a second incubation step, a monoclonal mouse-anti-PMN-Elastase antibody (AB) is added. This antibody is able to detect both the free and the complexed form with the specific inhibitor (α 1-Proteinase inhibitor = α 1-Antitrypsin). The quantification of the bound **PMN-Elastase** is carried out by adding an anti-mouse peroxidase-labeled conjugate (CONJ). Finally, the **PMN-Elastase-antigen-antibody-complex** is incubated with the peroxidase substrate, tetramethylbenzidine (SUB). An acidic stop solution (STOP) is then added to terminate the reaction. The color changes from blue to yellow. The intensity of the yellow color is directly proportional to the concentration of **PMN-Elastase** in the sample. A dose response curve of the absorbance unit (optical density, OD) vs. concentration is generated, using the values obtained from the calibrators. **PMN-Elastase**, present in the patient samples, is determined directly from this curve. The combination of two specific antibodies in the **PMN-Elastase** ELISA drastically reduces the possibility of false results and offers a reliable diagnostic system to the user.

Test procedure

1.	Prior to use in the assay allow all reagents and samples to come to room temperature and mix by gentle swirling and inversion.
2.	Mark the positions of STD (Standards)/ SAMPLE / CTRL (Control) in duplicate on a protocol sheet
3.	Take PLATE (microtiter strips) out of the kit. Store unused strips in the original package bag at 2-8° C. Strips are stable until expiry date stated on the label
4.	Wash the wells 5x with 250 µl of diluted wash buffer. After the last wash, remove remaining wash buffer by hitting the plate against paper towel
5.	Add 100 µl of STD (Standards)/ SAMPLE / CTRL (Control) in duplicate into respective well
6.	Cover the plate tightly and incubate for 1 hour at room temperature (18-26°C) on a horizontal mixer
7.	Aspirate and wash the wells 5x with 250 µl of diluted wash buffer. After the last wash, remove remaining wash buffer by hitting the plate against paper towel
8.	Add 100 µl of AB (Detection antibody, 2nd Antibody) into each wells mix gently

9.	Cover the plate tightly and incubate for 1 hour at room temperature (18-26°C) on a horizontal mixer
10.	Aspirate and wash the wells 5x with 250 µl of diluted wash buffer. After the last wash, remove remaining wash buffer by hitting the plate against paper towel
11.	Add 100 µl of CONJ (Conjugate) into each well
12.	Cover the plate tightly and incubate for 1 hour at room temperature (18-26°C) on a horizontal mixer
13.	Aspirate and wash the wells 5x with 250 µl of diluted wash buffer. After the last wash, remove remaining wash buffer by hitting the plate against paper towel
14.	Add 100 µl of SUB (Substrate) into each well
15.	Incubate for 5 - 20 minutes at room temperature (18-26°C) in the dark
16.	Add 50 µl of STOP (Stop solution) into each well, shake well
17.	Determine absorption immediately with an ELISA reader at 450 nm against 620 nm as a reference. If no reference wavelength is available, read only at 450 nm. If the extinction of the highest standard exceeds the range of the photometer, absorption must be measured immediately at 405 nm against 620 nm as a reference

8. RESULTS

The following algorithms can be used alternatively to calculate the results. We recommend to use the "4-Parameter-algorithm".

1. 4-parameter-algorithm

It is recommended to use a linear ordinate for optical density and a logarithmic abscissa for concentration. When using a logarithmic abscissa, the zero calibrator must be specified with a value less than 1 (e. g. 0.01).

2. Point-to-point-calculation

We recommend a linear ordinate for optical density and a linear abscissa for concentration.

3. Spline-algorithm

We recommend a linear ordinate for optical density and a logarithmic abscissa for concentration. When using a logarithmic abscissa, the zero calibrator must be specified with a value less than 1 (e. g. 0.01).

The plausibility of the pairs of values should be examined before the automatic evaluation of the results. If this option is not available with the used program, a control of the paired values should be done manually.

Faeces

To obtain the PMN-Elastase concentration in fecal samples multiply the estimated value by the dilution factor according to the sample preparation:

Sample preparation 1 and 2b: Dilution factor constant : 1 : 50

Multiply the obtained result by **50** to get the real concentration.

Sample preparation 2a: The dilution factor is variable,
dependent on the sample amount, and is taken from the table.

Multiply the obtained result by **the corresponding dilution factor** to get the real concentration.

Seminal plasma

For the calculation of the PMN-Elastase concentration in seminal plasma, the result has to be multiplied by **10 or 20** dependent on the sample dilution.

Serum

For the calculation of the PMN-Elastase concentration in serum, the result has to be multiplied by **500**.

Plasma

For the calculation of the PMN-Elastase concentration in plasma, the result has to be multiplied by **100**.

9. LIMITATIONS

Samples with PMN-Elastase levels greater than the highest calibrator, should be further diluted and re-assayed.

Increased concentrations of PMN-Elastase in seminal plasma support the diagnosis of an acute or chronic Adnexitis. In the case of simultaneous Urethritis anterior, the obtained results are elevated due to contamination of the ejaculate with granulocytes and can not be used. A possible presented Urethritis anterior is easy for diagnosis. In this case, the obtained results must be excluded or the patients treated before taking samples for PMN-Elastase analysis. Only such data can be used for a reliable diagnosis, especially by existing Adnexitis.

10. QUALITY CONTROL

Immundiagnostik recommends the use of commercial control samples for internal quality control.

Control samples or serum pools should be analyzed with each run of calibrators and patient samples. Results, generated from the analysis of control samples, should be evaluated for acceptability using appropriate statistical methods. The results for the patient sample may not be valid, if within the same assay one or more values of the quality control sample are lying outside the acceptable limits.

Expected values

(1g stool be equivalent to 1 ml)

PMN-Elastase concentrations in faeces of a healthy person (n=76):

< 62 ng/ml

PMN-Elastase concentrations in plasma of a healthy person (n=37):

19-78 ng/ml

PMN-Elastase concentrations in serum of a healthy person (n = 52):

average = 688 ng/ml (186 – 1991ng/ml)

We recommend each laboratory to establish its own norm concentration range.

11. PERFORMANCE CHARACTERISTICS

Precision and reproducibility

Intra-Assay (n=20)		
Probe	PMN Elastase [ng/ml]	Vk [%]
1	2.3	14
2	0.72	4.5

Inter-Assay (n=20)		
Probe	PMN Elastase [ng/ml]	Vk [%]
1	4.9	4
2	27.8	10

Recovery

Two samples were spiked with different PMN-Elastase calibrator amounts and measured with the assay.

Probe [ng/ml]	Spike [ng/ml]	PMN-Elastase Expected [ng/ml]	PMN-Elastase Measured [ng/ml]
1.5	3	4.5	4.9
	6	7.5	7.7
0.75	5	5.75	5.7
	1.65	2.4	2.7
	0.65	1.4	1.5

Sensitivity

The sensitivity limit was set as $B_0 + 2 \text{ SD}$. The Zero-standard was measured 20 times.

Sample	PMN-Elastase Mean value [OD]	Standard variation	Detection limit [ng/ml]
1	0.025	0.024	0.12

Cross reactivity

No cross reactivity to other proteins in stool and serum was observed.

Cross reactivity with PMN-elastase as well as a good correlation with PMN-elastase content in mouse serum was observed.

Linearity

Two patient samples were diluted with ELISA wash buffer and analyzed with the Immundiagnostik PMN-Elastase ELISA test.

Sample	Dilution	Expected [ng/ml]	Measured [ng/ml]
A	undiluted	22	22
	1:2	11	10.3
	1:4	5.5	5.3
	1:8	2.7	2.7
B	undiluted	4.6	4.6
	1:2	2.3	2.4
	1:4	1.2	1.3
	1:8	0.6	0.9

12. PRECAUTIONS

- For *in vitro* diagnostic use only.
- The quality control guidelines should be observed.
- Human material used in the kit components was tested and found to be negative for HIV, Hepatitis B and Hepatitis C. However, for safety reasons, all kit components should be treated as potentially infectious.
- Reagents of the kit package contain sodium azide and thimerosal as bactericides. Sodium azide and thimerosal are toxic. The substrates for the enzymatic color reactions are described to be also toxic and carcinogenic. Contact with skin or mucous membranes has to be avoided.
- Stop solution consists of sulfuric acid, a strong acid. Although diluted, it still must be handled with care. It can cause burns and should be handled with gloves, eye protection, and appropriate protective clothing. Any spill should be wiped out immediately with copious quantities of water.

13. TECHNICAL HINTS

- Do not interchange different lot numbers of any kit component within the same assay.
- Reagents should not be used beyond the expiration date shown on the kit label.
- Substrate solution should remain colourless until added to the plate.
- To ensure accurate results, proper adhesion of plate sealers during incubation steps is necessary.
- Avoid foaming when mixing reagents.
- The assay should always be performed according the enclosed manual.
- Incubation time, incubation temperature and the volumes of the different components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the test results. Immundiagnostik can therefore not be held responsible for any damage resulting from this.

14. GENERAL NOTES ON THE TEST AND TEST PROCEDURE

- This assay was produced and put on the market according to the IVD guidelines of 98/79/EC.
- The kit components made of human serum are tested for Hepatitis B, Hepatitis C and HIV and found to be negative. However, since no test method can offer complete assurance that infectious agents are absent, these reagents should be handled as recommended for any potentially infectious human serum or blood specimen. The normal precautions for laboratory work should be observed.
- All reagents in the kit package are for *in vitro* diagnostic use only.
- The guidelines for medical laboratories should be observed.
- Incubation time, incubation temperature and volumes of the different components are defined by the producer. Any variation of the test procedure, which is not coordinated with the producer, may influence the results of the test. Immundiagnostik can therefore not be held responsible for any damage resulting from this.
- Warranty claims and complaints in respect of deficiencies must be lodged within 14 days of received of the product. The product shall be send to Immundiagnostik together with the complaint in writing.

15. REFERENCES

1. Heinichen et al.: 1995; *Clin. Lab.* **41**, 539-545
2. Oremek et al.: 1995; *MTA* **10**, 273-278.

1/3/2008 1/3/2008_PMN-E.doc

Used symbols:

	Temperature limitation		Catalogue Number
	In Vitro Diagnostic Medical Device		Contains sufficient for <n> tests
	Manufacturer		Use by
	Lot number		