



# 암 치료의 새로운 가능성

R&D SYSTEMS  
a biotechne® brand

# Memory T Cells!



## MagCelect Human Memory CD4+ T Cell Isolation Kit (Cat# MAGH116)

- Magnetic Nanoparticle 기반 negative selection을 통해 간편 분리!
- Cell damage 없는 작은 크기의 beads 사용!
- 짧은 시간 내에 high purity 분리 가능!

구성품
Human CD4+ Memory T Cell Biotinylated Antibody Cocktail
Streptavidin Ferrofluid
10X Buffer

## Product List of Flow Cytometry Antibodies

Target	Species	APC	Fluorescein	PE	PerCP	488	700	405/594/64/7/750	Biotin	Unconjugated	Cat#
CD3ε	Hu, Mu	○	○	○		○	○	○/○/○/○		○	MAB100
CD28	Hu, Mu					○	○	○/○/○/○	○	○	MAB3422
CD3	Hu, Mu	○	○	○	○	○	○	○/○/○/○	○	○	MAB4841
CD4	Hu, Mu, Rt, Ca	○	○	○	○	○	○	○/○/○/○	○	○	MAB379
CD8α	Hu, Mu, Rt, Pr	○	○	○	○	○	○	○/○/○/○	○	○	MAB1509
CCR7	Hu, Mu	○	○	○		○	○	○/○/○/○	○	○	MAB197
CD27	Hu, Mu	○				○	○	○/○/○/○	○	○	MAB382
CD44	Hu, Mu, Rt, Ca	○		○		○	○	○/○/○/○		○	MAB6127
CD45RA	Hu					○	○	○/○/○/○		○	MAB11444
CD45RO	Hu			○		○	○	○/○/○/○		○	MAB10642
CD62L	Hu, Mu		○	○		○	○	○/○/○/○		○	MAB728
CD69	Hu, Mu	○		○		○	○	○/○/○/○	○	○	MAB23591
CD103	Mu					○			○	○	MAB1990
IFN-γ	Hu, Mu, Bv, Pr	○	○	○	○	○	○	○/○/○/○		○	MAB285
TNF-α	Hu, Mu		○	○		○	○	○/○/○/○	○	○	MAB9677
IL-2	Hu, Pr	○		○		○	○	○/○/○/○		○	MAB202
IL-10	Hu,					○	○	○/○/○/○		○	MAB217
B7-H1/PD-L1	Hu, Mu	○		○	○	○	○	○/○/○/○		○	MAB1561