



ELISA Kits for Cytokine Research

August 2010

- ✎ Chemokines & Receptors
- ✎ Common beta Chain Receptor Family
- ✎ Common gamma Chain Receptor Family
- ✎ EGF Family
- ✎ FGF Family
- ✎ IL-1 Ligands & IL-1R/TLR Family
- ✎ IL-6 Family
- ✎ IL-10/Interferon Family
- ✎ IL-17 Family
- ✎ Other Interleukins
- ✎ TGF-beta Superfamily
- ✎ TNF Superfamily

Chemokines & Receptors

Product	Species	Sample Types	Sensitivity	Range	Catalog #	Size
CCL2/JE/MCP-1	Canine	S P C	10.8 pg/mL	15.6-1,000 pg/mL	CACP00	1 Kit
	Human	S P C U	5 pg/mL	31.2-2,000 pg/mL	DCP00	1 Kit
	Mouse	S C	2 pg/mL	15.6-1,000 pg/mL	MJE00	1 Kit
CCL3/MIP-1 alpha	Human	S P C	10 pg/mL*	31.2-1,000 pg/mL*	DMA00	1 Kit
	Mouse	S C	1.5 pg/mL	4.7-300 pg/mL	MMA00	1 Kit
CCL4/MIP-1 beta	Human	S P C	11 pg/mL*	31.2-2,000 pg/mL*	DMB00	1 Kit
	Mouse	S P C	3 pg/mL	7.8-500 pg/mL	MMB00	1 Kit
CCL5/RANTES	Human	S P C U	6.6 pg/mL	31.2-2,000 pg/mL	DRN00B	1 Kit
	Mouse	S P C	2 pg/mL	7.8-500 pg/mL	MMR00	1 Kit
CCL7/MCP-3/MARC	Human	S P C	8.52 pg/mL	15.6-1,000 pg/mL	DCC700	1 Kit
CCL11/Eotaxin	Human	S P C	5 pg/mL	15.6-1,000 pg/mL	DTX00	1 Kit
	Mouse	S P C	3 pg/mL	15.6-1,000 pg/mL	MME00	1 Kit
CCL12/MCP-5	Mouse	S P C	2.1 pg/mL	15.6-1,000 pg/mL	MCC120	1 Kit
CCL17/TARC	Human	S P C	7 pg/mL	31.2-2,000 pg/mL	DDN00	1 Kit
	Mouse	S P C	5 pg/mL	7.8-500 pg/mL	MCC170	1 Kit
CCL20/MIP-3 alpha	Human	S P C	0.87 pg/mL	7.8-500 pg/mL	DM3A00	1 Kit
	Mouse	S P C	0.49 pg/mL	3.9-250 pg/mL	MCC200	1 Kit
CCL21/6Ckine	Human	S P C	33.5 pg/mL	78.1-5,000 pg/mL	D6C00	1 Kit
CCL22/MDC	Human	S P C	62.5 pg/mL	125-4,000 pg/mL	DMD00	1 Kit
	Mouse	S P C	1.8 pg/mL	7.8-500 pg/mL	MCC220	1 Kit
CCL24/Eotaxin-2/MPIF-2	Human	S P C	14.3 pg/mL	78-5,000 pg/mL	DCC240B	1 Kit
CCL26/Eotaxin-3	Human	S P C	5.2 pg/mL	7.8-500 pg/mL	DCC260	1 Kit
CCL27/CTACK	Human	S P C	4.68 pg/mL	15.6-1,000 pg/mL	DCC270	1 Kit
Chemerin	Mouse	S P C	3.47 pg/mL	46.9-3,000 pg/mL	MCHM00	1 Kit
CX3CL1/Fractalkine	Human	S P C U V M	0.072 ng/mL	0.156-10 ng/mL	DCX310	1 Kit
	Mouse	S P C T	0.32 ng/mL	0.62-40 ng/mL	MCX310	1 Kit
CXCL1/GRO alpha/KC/CINC-1	Human	S P C	10 pg/mL	31.2-1,000 pg/mL	DGR00	1 Kit
	Mouse	S C	2 pg/mL	15.6-1,000 pg/mL	MKC00B	1 Kit
	Rat	S P C	1.3 pg/mL	7.8-500 pg/mL	RCN100	1 Kit
CXCL2/GRO beta/MIP-2/CINC-3	Mouse	S C	1.5 pg/mL	7.8-500 pg/mL	MM200	1 Kit
	Rat	S P C	2.7 pg/mL	31.2-2,000 pg/mL	RCN300	1 Kit
CXCL3/GRO gamma/CINC-2/DCIP-1	Rat	S P C	3.9 pg/mL	6.25-400 pg/mL	RCN200	1 Kit
CXCL5/ENA-78	Human	S P C	15 pg/mL	31.2-2,000 pg/mL	DX000	1 Kit
CXCL6/GCP-2	Human	S P C	8 pg/mL	31.2-2,000 pg/mL	DGC00	1 Kit
CXCL8/IL-8	Canine	S P C	4.31 pg/mL	15.6-1,000 pg/mL	CA8000	1 Kit
	Human	S P C	7.5 pg/mL	31.2-2,000 pg/mL	D8000C	1 Kit
	Porcine	S C	6.7 pg/mL	62.5-4,000 pg/mL	P8000	1 Kit
CXCL8/IL-8 Chemiluminescent	Human	S P C	0.97 pg/mL	1.6-5,000 pg/mL	Q8000B	1 Kit
CXCL9/MIG	Human	S P C	11.3 pg/mL	31.2-2,000 pg/mL	DCX900	1 Kit
	Mouse	S P C T	7.8 pg/mL	31.2-2,000 pg/mL	MCX900	1 Kit
CXCL10/IP-10/CRG-2	Human	S P C V	4.46 pg/mL	7.8-500 pg/mL	DIP100	1 Kit
	Mouse	S P C	4.2 pg/mL	31.2-2,000 pg/mL	MCX100	1 Kit
CXCL11/I-TAC	Human	S P C	39.7 pg/mL	62.5-4,000 pg/mL	DCX110	1 Kit
CXCL12/SDF-1 alpha	Human	P C	47 pg/mL	156-10,000 pg/mL	DSA00	1 Kit

For research use only. Not for use in diagnostic procedures unless noted. Economy packs available, please inquire.

*Sample-specific, please refer to product package insert for sensitivities and ranges for all sample types.

Sample Types:

S - Serum C - Cell Culture Supernate V - Saliva E - Cell Extracts L - Cell Lysate PL - Platelet Lysate T - Tissue Homogenate
P - Plasm U - Urine B - Bone Extracts F - Follicular Fluid M - Breast Milk SF - Subcellular Fractions

Chemokines & Receptors

Product	Species	Sample Types	Sensitivity	Range	Catalog #	Size
CXCL12/SDF-1 alpha	Mouse	S C	0.069 ng/mL	0.156-10 ng/mL	MCX120	1 Kit
CXCL13/BLC/BCA-1	Human	S P C V	3.97 pg/mL	7.8-500 pg/mL	DCX130	1 Kit
	Mouse	S P C	2.84 pg/mL	15.6-1,000 pg/mL	MCX130	1 Kit
CXCL16	Human	S P C	0.017 ng/mL	0.156-10 ng/mL	DCX160	1 Kit

Common beta Chain Receptor Family

Product	Species	Sample Types	Sensitivity	Range	Catalog #	Size
GM-CSF	Human	S P C	3 pg/mL	7.8-500 pg/mL	DGM00	1 Kit
	Mouse	S C	5.8 pg/mL	7.8-500 pg/mL	MGM00	1 Kit
	Rat	S C	2.1 pg/mL	15.6-1,000 pg/mL	RGM00	1 Kit
GM-CSF High Sensitivity	Human	S P U	0.26 pg/mL*	0.5-32 pg/mL*	HSGM0	1 Kit
IL-3	Human	S P C	7.4 pg/mL	31.2-2,000 pg/mL	D3000	1 Kit
	Mouse	S C	2.5 pg/mL	7.8-500 pg/mL	M3000	1 Kit
IL-5	Human	S P C U	1.08 pg/mL	3.9-250 pg/mL	D5000B	1 Kit
	Mouse	S C	7 pg/mL	15.6-1,000 pg/mL	M5000	1 Kit

Common gamma Chain Receptor Family

Product	Species	Sample Types	Sensitivity	Range	Catalog #	Size
IL-2	Human	S P C	7 pg/mL	31.2-2,000 pg/mL	D2050	1 Kit
	Mouse	S C	3 pg/mL	15.6-1,000 pg/mL	M2000	1 Kit
	Rat	S P C	15 pg/mL	31.2-2,000 pg/mL	R2000	1 Kit
IL-2 Chemiluminescent	Human	S P C	0.25 pg/mL	1.7-1,250 pg/mL	Q2000B	1 Kit
IL-2 R alpha	Human	S P C	10 pg/mL	78-5,000 pg/mL	DR2A00	1 Kit
IL-4	Human	S P C	10 pg/mL	31.2-2,000 pg/mL	D4050	1 Kit
	Mouse	S C	2 pg/mL	7.8-500 pg/mL	M4000B	1 Kit
	Rat	S C	5 pg/mL	15.6-1,000 pg/mL	R4000	1 Kit
IL-4 Chemiluminescent	Human	S P C	1 pg/mL	1.6-5,000 pg/mL	Q4000	1 Kit
IL-4 High Sensitivity	Human	S C	0.22 pg/mL	0.25-16 pg/mL	HS400	1 Kit
IL-7	Mouse	S P C	8.3 pg/mL	31.2-2,000 pg/mL	M7000	1 Kit
IL-7 High Sensitivity	Human	S P C	0.1 pg/mL*	0.156-10 pg/mL*	HS750	1 Kit
IL-15	Human	S P C	2 pg/mL	3.9-250 pg/mL	D1500	1 Kit
	Human	S P C	0.431 pg/mL	1.03-750 pg/mL	Q1500B	1 Kit

EGF Family

Product	Species	Sample Types	Sensitivity	Range	Catalog #	Size
EGF	Human	S P C U	0.7 pg/mL	3.9-250 pg/mL	DEG00	1 Kit
	Mouse	S P C U T	1.64 pg/mL	7.8-500 pg/mL	MEG00	1 Kit
EGF R/ErbB1	Human	S P C M	0.036 ng/mL	0.313-20 ng/mL	DEGFR0	1 Kit
Phospho-EGF R/ErbB1 (Y1068), Cell-Based ELISA	Human	WC			KCB1095	1 Kit

For research use only. Not for use in diagnostic procedures unless noted. Economy packs available, please inquire.

*Sample-specific, please refer to product package insert for sensitivities and ranges for all sample types.

Sample Types:

S - Serum C - Cell Culture Supernate V - Saliva E - Cell Extracts L - Cell Lysate PL - Platelet Lysate T - Tissue Homogenate
P - Plasm U - Urine B - Bone Extracts F - Follicular Fluid M - Breast Milk SF - Subcellular Fractions

EGF Family

Product	Species	Sample Types	Sensitivity	Range	Catalog #	Size
Phospho-ErbB2/Her2 (Y1196), Cell-Based ELISA	Human	WC			KCB4438	1 Kit
Phospho-ErbB3/Her3 (Y1262), Cell-Based ELISA	Human	WC			KCB5677	1 Kit
TGF-alpha	Human	S P C M	7.1 pg/mL	15.6-1,000 pg/mL	DTGA00	1 Kit

FGF Family

Product	Species	Sample Types	Sensitivity	Range	Catalog #	Size
FGF acidic	Human	S P C U	13.9 pg/mL	31.2-2,000 pg/mL	DFA00B	1 Kit
FGF basic	Human	S P C	3 pg/mL	10-640 pg/mL	DFB50	1 Kit
	Human	S P U	0.091 pg/mL	0.313-20 pg/mL	HSFB00D	1 Kit
KGF/FGF-7	Human	S P C	15 pg/mL	31.2-2,000 pg/mL	DKG00	1 Kit
FGF-19	Human	S P C	3.35 pg/mL	15.6-1,000 pg/mL	DF1900	1 Kit
FGF-21	NEW Human	S P C	8.69 pg/mL	31.3-2,000 pg/mL	DF2100	1 Kit
	Mouse	S P C	13.4 pg/mL	31.3-2,000 pg/mL	MF2100	1 Kit
Phospho-FRS2 (Y436), Cell-Based ELISA	Human/Mouse/Rat	WC			KCB5126	1 Kit
Pentraxin 3/TSG-14	Human	P C V	0.116 ng/mL	0.31-20 ng/mL	DPTX30	1 Kit
	Mouse	S P C	0.02 ng/mL	0.23-15 ng/mL	MPTX30	1 Kit

IL-1 Ligands & IL-1R/TLR Family

Product	Species	Sample Types	Sensitivity	Range	Catalog #	Size
IL-1 alpha/IL-1F1	Human	S P C	1 pg/mL	3.9-250 pg/mL	DLA50	1 Kit
	Mouse	S C	2.5 pg/mL	4.69-300 pg/mL	MLA00	1 Kit
	Rat	S P C	4.12 pg/mL	15.6-1,000 pg/mL	RRA00	1 Kit
IL-1 beta/IL-1F2	Human	S P C	1 pg/mL	3.9-250 pg/mL	DLB50	1 Kit
	Mouse	S C	3 pg/mL	7.8-500 pg/mL	MLB00B	1 Kit
	Porcine	C	10 pg/mL*	19.5-1,250 pg/mL*	PLB00	1 Kit
	Rat	S P C	5 pg/mL	31.2-2,000 pg/mL	RLB00	1 Kit
IL-1 beta/IL-1F2 Chemiluminescent	Human	S P C	0.55 pg/mL	1.4-1,000 pg/mL	QLB00B	1 Kit
IL-1 beta/IL-1F2 High Sensitivity	Human	S P	0.14 pg/mL	0.125-8 pg/mL	HSLB00C	1 Kit
Pro IL-1 beta/IL-1F2	Human	S P C	7.2 pg/mL*	23.4-1,500 pg/mL*	DLBP00	1 Kit
IL-1ra/IL-1F3	Human	S P C	18.3 pg/mL	31.2-2,000 pg/mL	DRA00B	1 Kit
	Mouse	S P C	13 pg/mL	31.2-2,000 pg/mL	MRA00	1 Kit
IL-1 RII	Human	S P C	10 pg/mL	31.2-2,000 pg/mL	DR1B00	1 Kit
IL-18/IL-1F4	Human	S P	12.5 pg/mL	25.6-1,000 pg/mL	7620	1 Kit
	Mouse	S	25 pg/mL	25.6-1,000 pg/mL	7625	1 Kit
IL-33	Mouse	S P C	14.3 pg/mL	31.3-2,000 pg/mL	M3300	1 Kit
ST2/IL-1 R4	Human	S P C	13.5 pg/mL	31.3-2,000 pg/mL	DST200	1 Kit

For research use only. Not for use in diagnostic procedures unless noted. Economy packs available, please inquire.

*Sample-specific, please refer to product package insert for sensitivities and ranges for all sample types.

Sample Types:

S - Serum C - Cell Culture Supernate V - Saliva E - Cell Extracts L - Cell Lysate PL - Platelet Lysate T - Tissue Homogenate
 P - Plasm U - Urine B - Bone Extracts F - Follicular Fluid M - Breast Milk SF - Subcellular Fractions

IL-6 Family

Product	Species	Sample Types	Sensitivity	Range	Catalog #	Size
CNTF	Human	S P C	8 pg/mL	31.2-2,000 pg/mL	DNT00	1 Kit
G-CSF	Human	S P C	20 pg/mL	39-2,500 pg/mL	DCS50	1 Kit
	NEW Human	S P C	4.62 pg/mL	4.69-300 pg/mL	HSTCS0	1 Kit
	Mouse	S C	5 pg/mL	14.1-900 pg/mL	MCS00	1 Kit
gp130	Human	S P C	0.05 ng/mL*	0.125-8 ng/mL*	DGP00	1 Kit
IL-6	Canine	S P C	11.8 pg/mL	31.3-2,000 pg/mL	CA6000	1 Kit
	Human	S P C	0.7 pg/mL	3.12-300 pg/mL	D6050	1 Kit
	Mouse	S P C	1.8 pg/mL	7.8-500 pg/mL	M6000B	1 Kit
	Porcine	S P C	10 pg/mL	39.1-2,500 pg/mL	P6000	1 Kit
	Rat	S P C	36 pg/mL	62.5-4,000 pg/mL	R6000B	1 Kit
IL-6 Chemiluminescent	Human	S P C	0.35 pg/mL	0.48-1,500 pg/mL	Q6000B	1 Kit
IL-6 High Sensitivity	Human	S P U	0.11 pg/mL	0.156-10 pg/mL	HS600B	1 Kit
IL-6 R alpha	Human	S P C U	15.1 pg/mL	31.2-2,000 pg/mL	DR600	1 Kit
IL-11	Human	S P C	8 pg/mL*	15.6-1,000 pg/mL*	D1100	1 Kit
	Human	S P C	7.8 pg/mL	15.6-1,000 pg/mL	DLP00	1 Kit
Leptin/OB	Human	S P C	7.8 pg/mL	15.6-1,000 pg/mL	DLP00	1 Kit
	Mouse	S P C	22 pg/mL	62.5-4,000 pg/mL	MOB00	1 Kit
Leptin R	Human	S P C	0.128 ng/mL	0.312-20 ng/mL	DOBR00	1 Kit
LIF	Human	S P C	8 pg/mL	31.2-2,000 pg/mL	DLF00	1 Kit
	Mouse	S P C	3.13 pg/mL	21.9-1,400 pg/mL	MLF00	1 Kit
Phospho-STAT3 (Y705), Cell-Based ELISA	Human/Mouse	WC			KCB4607	1 Kit

IL-10/Interferon Family

Product	Species	Sample Types	Sensitivity	Range	Catalog #	Size
IFN-alpha	Human	C	12.5 pg/mL	10-500 pg/mL	41100-1	1 Kit
	e Date Mouse				42120-1-1	1 Kit
IFN-alpha Multi-Subtype	Human	C	12.5 pg/mL	10-500 pg/mL	41105-1	1 Kit
IFN-alpha Serum Sample	Human	S C	12.5 pg/mL	10-500 pg/mL	41110-1	1 Kit
IFN-beta	Human	C	25 pg/mL	25-2,000 pg/mL	41410-1	1 Kit
	Mouse	S C	15.5 pg/mL	15.6-1,000 pg/mL	42400-1	1 Kit
IFN-gamma	Canine	S P C	60 pg/mL	62.5-4,000 pg/mL	CAIF00	1 Kit
	Human	S P C	8 pg/mL	15.6-1,000 pg/mL	DIF50	1 Kit
	Mouse	S C	2 pg/mL	9.4-600 pg/mL	MIF00	1 Kit
	Porcine	S P C	11.2 pg/mL	39-2,500 pg/mL	PIF00	1 Kit
	Rat	S C	10 pg/mL	31.2-2,000 pg/mL	RIF00	1 Kit
IFN-omega	Human	S P C	5 pg/mL	3-500 pg/mL	41395-1	1 Kit
IL-10	Canine	S P C	3.8 pg/mL	15.6-1,000 pg/mL	CA1000	1 Kit
	Human	S P C	3.9 pg/mL	7.8-500 pg/mL	D1000B	1 Kit
	e Date Human	S P	0.45 pg/mL	0.78-50 pg/mL	HS100C	1 Kit
	Mouse	S C	4 pg/mL	15.6-1,000 pg/mL	M1000	1 Kit
	Porcine	S P C	5.5 pg/mL	31.2-2,000 pg/mL	P1000	1 Kit
	Rat	S P C	10 pg/mL	31.2-2,000 pg/mL	R1000	1 Kit
IL-10 Chemiluminescent	Human	S P C	3.03 pg/mL	3.2-10,000 pg/mL	Q1000	1 Kit
IL-19	Human	S P C U M	12.2 pg/mL	31.2-2,000 pg/mL	D1900	1 Kit

For research use only. Not for use in diagnostic procedures unless noted. Economy packs available, please inquire.

*Sample-specific, please refer to product package insert for sensitivities and ranges for all sample types.

Sample Types:

S - Serum C - Cell Culture Supernate V - Saliva E - Cell Extracts L - Cell Lysate PL - Platelet Lysate T - Tissue Homogenate
P - Plasm U - Urine B - Bone Extracts F - Follicular Fluid M - Breast Milk SF - Subcellular Fractions

IL-10/Interferon Family

Product	Species	Sample Types	Sensitivity	Range	Catalog #	Size
IL-20	Human	S P C	16.6 pg/mL	62.5-4,000 pg/mL	DL200	1 Kit
	Mouse	S P C	6.4 pg/mL	15.6-1,000 pg/mL	ML200	1 Kit
IL-22	Human	S P C U	5.8 pg/mL	15.6-1,000 pg/mL	D2200	1 Kit
	Mouse/Rat	S P C	8.2 pg/mL	15.6-1,000 pg/mL	M2200	1 Kit
Phospho-STAT2 (Y689), Cell-Based ELISA	Human	WC			KCB2890	1 Kit
Phospho-STAT3 (Y705), Cell-Based ELISA	Human/Mouse	WC			KCB4607	1 Kit
Phospho-STAT5a/b (Y699), Cell-Based ELISA	Human/Mouse	WC			KCB4190	1 Kit

IL-17 Family

Product	Species	Sample Types	Sensitivity	Range	Catalog #	Size
IL-17/IL-17A	Human	S P C	15 pg/mL	31.2-2,000 pg/mL	D1700	1 Kit
	Mouse	S P C	5 pg/mL	10.9-700 pg/mL	M1700	1 Kit
IL-17A/F Heterodimer	Mouse	S P C	2.8 pg/mL	7.81-500 pg/mL	M17AF0	1 Kit

Other Interleukins

Product	Species	Sample Types	Sensitivity	Range	Catalog #	Size
IL-12	Human	S P C U	5 pg/mL	7.8-500 pg/mL	D1200	1 Kit
IL-12 High Sensitivity	Human	S P C	0.5 pg/mL*	0.625-40 pg/mL*	HS120	1 Kit
IL-12 p70	Mouse	S C	2.5 pg/mL	7.8-500 pg/mL	M1270	1 Kit
IL-12/IL-23 p40	Porcine	S P C	18.2 pg/mL	46.9-3,000 pg/mL	P1240	1 Kit
IL-12/IL-23 p40 Allele-specific	Mouse	S C	4 pg/mL	11.7-750 pg/mL	M1240	1 Kit
IL-12/IL-23 p40 Homodimer	Human	S P C	15 pg/mL	31.2-2,000 pg/mL	DP400	1 Kit
IL-12/IL-23 p40 Nonallele-specific	Mouse	S P C	2.7 pg/mL	11.7-750 pg/mL	MP400	1 Kit
IL-13	Human	S P C	32 pg/mL	62.5-4,000 pg/mL	D1300	1 Kit
	Mouse	C	1.5 pg/mL	7.8-500 pg/mL	M1300CB	1 Kit
IL-16	Human	S P C	13.4 pg/mL	31.2-2,000 pg/mL	D1600	1 Kit
IL-23	Human	S P C	16.3 pg/mL	39-2,500 pg/mL	D2300B	1 Kit
	Mouse	S P C	4.17 pg/mL	15.6-1,000 pg/mL	M2300	1 Kit
IL-27 p28	Mouse	S P C	4.7 pg/mL	7.8-500 pg/mL	M2728	1 Kit

TGF-beta Superfamily

Product	Species	Sample Types	Sensitivity	Range	Catalog #	Size
Activin A	Human/Mouse/Rat	S P C V	7.85 pg/mL	15.6-1,000 pg/mL	DAC00B	1 Kit
BMP-2	Human/Mouse/Rat	C B	29 pg/mL	62.5-4,000 pg/mL	DBP200	1 Kit
BMP-4	Human	C B	3.68 pg/mL	31.2-2,000 pg/mL	DBP400	1 Kit
BMP-7	Human	S P C U B	7.83 pg/mL	31.2-2,000 pg/mL	DBP700	1 Kit
Endoglin/CD105	Human	S P C	0.03 ng/mL	0.156-10 ng/mL	DNDG00	1 Kit

For research use only. Not for use in diagnostic procedures unless noted. Economy packs available, please inquire.

*Sample-specific, please refer to product package insert for sensitivities and ranges for all sample types.

Sample Types:

S - Serum C - Cell Culture Supernate V - Saliva E - Cell Extracts L - Cell Lysate PL - Platelet Lysate T - Tissue Homogenate
P - Plasm U - Urine B - Bone Extracts F - Follicular Fluid M - Breast Milk SF - Subcellular Fractions

TGF-beta Superfamily

Product	Species	Sample Types	Sensitivity	Range	Catalog #	Size
Follistatin	Human	S P C F	83 pg/mL	250-16,000 pg/mL	DFN00	1 Kit
GDF-15	Human	S P C U	4.39 pg/mL	23.4-1,500 pg/mL	DGD150	1 Kit
TGF-beta 1	Human	S P C U	15.4 pg/mL	31.2-2,000 pg/mL	DB100B	1 Kit
	Mouse/Rat/Porcine/ Canine	S P C	15.4 pg/mL	31.2-2,000 pg/mL	MB100B	1 Kit
TGF-beta 2	Human	S P C	7 pg/mL	31.2-2,000 pg/mL	DB250	1 Kit

TNF Superfamily

Product	Species	Sample Types	Sensitivity	Range	Catalog #	Size
BAFF/BlyS/TNFSF13B	Human	S P C	11.9 pg/mL*	62.5-4,000 pg/mL*	DBLYS0	1 Kit
	Mouse	S P C	7.8 pg/mL	46.9-3,000 pg/mL	MBLYS0	1 Kit
CD40 Ligand/TNFSF5	Human	S P C	10.1 pg/mL	62.5-4,000 pg/mL	DCDL40	1 Kit
Fas/TNFRSF6/CD95	Human	S P C	20 pg/mL	31.2-2,000 pg/mL	DFS00	1 Kit
Fas Ligand/TNFSF6	Human	S P C L	8.05 pg/mL	15.6-1,000 pg/mL	DFL00	1 Kit
	Mouse	S P C	3.9 pg/mL	31.2-2,000 pg/mL	MFL00	1 Kit
LIGHT/TNFSF14	Human	S P C V	16.5 pg/mL	31.2-2,000 pg/mL	DLIT00	1 Kit
Osteoprotegerin/TNFRSF11B	Mouse	S P C	4.5 pg/mL	31.2-2,000 pg/mL	MOP00	1 Kit
TNF-alpha/TNFSF1A	Canine	S P C	4.2 pg/mL	7.8-500 pg/mL	CATA00	1 Kit
	Human	S P C	5.5 pg/mL	15.6-1,000 pg/mL	DTA00C	1 Kit
	Mouse	S C	5.1 pg/mL	23.4-1,500 pg/mL	MTA00	1 Kit
	Porcine	S P C	5 pg/mL	23.4-1,500 pg/mL	PTA00	1 Kit
	Rat	S P C	5 pg/mL	12.5-800 pg/mL	RTA00	1 Kit
	Rhesus Macaque	S P C	6.2 pg/mL	15.6-1,000 pg/mL	RHMTA0	1 Kit
TNF-alpha/TNFSF1A Chemiluminescent	Human	S P C	0.736 pg/mL	2.2-7,000 pg/mL	QTA00B	1 Kit
TNF-alpha/TNFSF1A High Sensitivity	Human	S P	0.191 pg/mL	0.5-32 pg/mL	HSTA00D	1 Kit
TNF RI/TNFRSF1A	Human	S P C U	1.2 pg/mL	7.8-500 pg/mL	DRT100	1 Kit
	Mouse	S C	5 pg/mL	15.6-1,000 pg/mL	MRT10	1 Kit
TNF RII/TNFRSF1B	Human	S P C U	2.3 pg/mL	7.8-500 pg/mL	DRT200	1 Kit
	Mouse	S P C	5 pg/mL	15.6-1,000 pg/mL	MRT20	1 Kit
TRAIL/TNFSF10	Human	S P C V L	7.87 pg/mL	15.6-1,000 pg/mL	DTRL00	1 Kit
TRANCE/TNFSF11/RANK L	Mouse	S C	5 pg/mL	31.2-2,000 pg/mL	MTR00	1 Kit

For research use only. Not for use in diagnostic procedures unless noted. Economy packs available, please inquire.

*Sample-specific, please refer to product package insert for sensitivities and ranges for all sample types.

Sample Types:

S - Serum C - Cell Culture Supernate V - Saliva E - Cell Extracts L - Cell Lysate PL - Platelet Lysate T - Tissue Homogenate
P - Plasm U - Urine B - Bone Extracts F - Follicular Fluid M - Breast Milk SF - Subcellular Fractions