

Biosimilar activity test를 위한 tool!

Biosimilar Reference

Biological Characterization of SB3, a Trastuzumab Biosimilar, and the Influence of Changes in Reference Product Characteristics on the Similarity Assessment

[Click here >>](#)


Demonstration of functional similarity of a biosimilar adalimumab SB5 to Humira®

[Click here >>](#)

In-Use Stability of the Rituximab Biosimilar CT-P10 (Truxima®) Following Preparation for Intravenous Infusion and Storage

[Click here >>](#)

Quidel C1q Protein for C1q Binding Test



Biosimilar가 효과적으로 cancer cell death를 유도하기 위해서는 C1q와의 binding 효율이 중요합니다. Quidel에서는 functional activity가 확인된 human C1q protein을 제공함으로써 제작된 biosimilar의 C1q binding 효율을 효과적으로 확인할 수 있습니다.

Order information


Cat#	Name
A400	C1q Protein (1 ML)

Quidel Normal Human Serum Complement for CDC Test

Biosimilar가 실제 cancer cell death를 얼마나 유도하는지 in vitro에서 확인하기 위해서는 체내에 존재하는 complement를 첨가해주어야 합니다. Quidel에서는 CDC test에 사용할 수 있는 normal human serum complement를 제공하고 있습니다. 제품 출시 전에 면역반응 여부를 체크함으로써 효과적인 CDC test가 이루어 질 수 있도록 합니다.

Order information

Cat#	Name
A112	Normal Human Serum Complement (2.5 ML)
A113	Normal Human Serum Complement (5 ML)



- Normal Human Serum Complement
- Biosimilar
- Cancer cell